

Quick Start Guide

CORTEXA Note | Enjoy Caring Again



Scan to download the Cortexa App.
Activation is required before first use.

1. Power On & Pair


Press and hold the power button. The LED will blink white, indicating the device is in pairing mode.

Download the Cortexa app from the App Store or Google Play. Create your account and follow the instructions to complete pairing. Once connected, the LED will turn solid white.



2. Start Recording

Clip the device to your coat. Press the record button once to start. The LED will turn solid red, confirming recording is active.



3. End & Transfer

Press the record button again to stop. The audio will automatically sync to your mobile app via Bluetooth or Wi-Fi. Please keep the app open during sync.



4. Review & Integrate

Review your AI-generated SOAP notes and Patient Instructions instantly on the mobile app or at app.cortexanote.com.



- **Clinical-Grade Hardware**

Purpose-built wearable with dual-mic noise cancellation, ensuring clear capture even in busy hospital environments.



- **Cross-Platform Access**

Access your clinical notes anywhere through our iOS, Android, and web applications with seamless sync.



- **HIPAA-Compliant Security**

HIPAA-compliant platform featuring end-to-end encryption and secure cloud processing.



- **Seamless EHR Integration**

One-click plugin import for supported systems, or universal copy-paste for any EHR.



- **Medical-Specific AI**

Powered by specialized medical AI. Trained on millions of medical records to accurately capture complex terminology.



- **Smart Medical Library**

Tailored for 144+ Specialties. Access precision terminology and custom templates designed for your specific field.



FCC compliance statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Radio Frequency (RF) Exposure Information

This device is designed to meet the requirements for exposure to radio waves established by the FCC. The device complies with FCC RF exposure limits for general population/uncontrolled exposure. To maintain compliance, please use the device as described in the user manual.

Warranty Policy

Your Cortexa Note comes with a 12-month worry-free warranty (or as required by local laws) starting from the date of purchase. Cortexa AI, Inc. will provide free repair or replacement for performance failures under normal use. Please retain your proof of purchase.

This warranty does not cover:

- Unauthorized disassembly, repair, or modification.
- Damage resulting from misuse, drops, liquid exposure, or improper accessory use.
- Damage caused by force majeure (e.g., natural disasters).
- Normal wear and tear or cosmetic damage.
- Products with altered/missing serial numbers or those purchased from unauthorized resellers.

For support, contact us at: support@cortexanote.com

Clinical Disclaimer: Cortexa Note is an assistive documentation tool. The clinician is responsible for reviewing and verifying the accuracy of all generated notes and medical orders before entry into the EHR. Cortexa AI, Inc. is not responsible for medical errors.

Download the Cortexa Note app



Visit cortexanote.com to unlock full features

We're Here to Help

Have questions?

Our support team is ready to assist.

Visit support.cortexanote.com or scan the QR code for FAQs and troubleshooting guides.

